Fibers & Probes



| 138 | Overview |
|-----|---|
| 139 | Custom Fiber & Probe Options |
| 142 | Premium-grade Assemblies |
| 143 | Unjacketed Bulk Optical Fiber |
| 144 | Laboratory-grade Assemblies |
| 146 | Xtreme Solarization-resistant Assemblies |
| 147 | Flame Loop Fiber Optic Probe |
| 147 | Cosine-corrected Irradiance Probe |
| 148 | Reflection/Backscattering Probes |
| 150 | Transmission Dip Probes |
| 152 | Industrial Process Probe |
| 153 | Vacuum Feedthroughs |
| 154 | Optical Fiber Kits |
| 157 | Fiber & Probe Fixtures & Holders |
| 158 | Fiber & Probe Accessories |

Overview: Fibers & Probes

The Most Flexible Line in the Industry


Ocean Optics is the most versatile supplier of optical fibers and accessories for spectroscopy in the industry. We offer everything from one-off patch cords and custom assemblies, to OEM builds for all sorts of applications. Our fiber accessories, fixtures and fiber assembly kits allow you to easily connect or manipulate fibers, and integrate them into tricky experiment set-ups. Optical fiber technology has been paramount to our success and makes possible our "take the instrument to the sample" maxim. Fiber is the nucleus of our analytical instrumentation and accessory design philosophy. And it's optical fiber technology that helped us create the world's most flexible line of sensing instruments.

Anatomy of an Assembly

At the fiber's core is pure silica; it's the diameter of the core that you need to consider when purchasing an optical fiber assembly. (The core diameter is often in the product's item code. For example, the P600-UV-VIS has a 600 μ m diameter silica core.) Surrounding the core is a doped-fluorine silica cladding. A buffer material is then applied. A buffer coats the core and cladding , strengthens the fiber and reduces stray light even further. In most assemblies polyimide is used as the buffer; other assemblies use aluminum or acrylate. Then a jacketing is applied over the core, cladding and buffer to protect the fiber and provide strain relief. For off-the-shelf Premium-grade "Q" Optical Fiber Assemblies, the standard jacketing is stainless steel silicone monocoil. For off-the-shelf the Laboratory-grade Optical Fiber Assemblies, the standard jacketing is zip tube blue PVDF. (There are several other jacketing options when creating a custom assembly.) Precision SMA 905 Connectors terminate the assembly and are precisely aligned to the spectrometer's slit to ensure concentricity of the fiber. Finally, captive end caps protect the fiber tips against scratches and contaminants.

Assembly Identifiers

Our optical fiber and probe assemblies are clearly and cleanly labeled in three ways so that you always know the following about your assembly: its name, its core diameter, and its most efficient wavelength region.



BOOT COLLARS

The assembly's boot color lets you know the fiber type and the most efficient wavelength range in which your fiber will work.

BAND COLORS

A color band tells you the diameter fiber with which you are working.




| Boot Color | Fiber Type | Most Efficient Wavelength Range | Premuim-grade Optical Fiber Assembly for each Fiber Type |
|---------------|-----------------------------------|------------------------------------|---|
| Gray | UV-VIS XSR Solarization-resistant | 180-900 nm | oceanoptics.com |
| Gray | UV/SR-VIS High OH content | 200-1100 nm | oceanoptics.com |
| Blue | UV-VIS High OH content | 300-1100 nm | oceanoptics.com |
| Red | VIS-NIR Low OH content | 400-2500 nm | oceanoptics.com |

Custom Fiber & Probe Assemblies

Custom Fiber & Probe Quoting Website


Continuous innovation is a critical ingredient for our growth. Our parent company, Halma p.l.c., believes that innovation is not just the responsibility of our development departments but is integral to all commercial activities within the business. Halma's successful Innovation Initiative allows all employees to deliver innovative ideas to help their companies achieve growth objectives. In 2006, Halma awarded Ocean Optics employees for our Custom Fiber Assembly Quoting and CAD Design website, designed to rapidly create quotations and CAD drawings of custom optical fiber assemblies.


The website allows our application sales scientists and customer service and support personnel to custom-create an optical fiber or probe assembly that automatically generates quotes and drawings for customers to view, modify and approve. Customers have hundreds of options available to them when creating a custom assembly. Our team of engineers uses our award-winning site to help guide you through choosing the best options for your unique application.





Custom Option: Fiber Type

One important consideration in building a custom assembly is which fiber type you should specify for your application. Typically, the wavelength range needed for your application should match the wavelength range of the fiber type. Here are the attenuation spectra for each fiber type we offer.











Custom Fiber & Probe Assemblies

Custom Option: Optical Fiber Diameter Sizes

After selecting the best Fiber Type, you should consider the diameter size of the pure silica core needed inside of your assembly. We offer several diameter sizes, and can recommend the appropriate assembly based on these criteria:

- 1. How much light do you need for your application? Reflection and fluorescence applications generally need more light, and larger diameter fibers are often better choices than smaller diameter fibers. For a laser application, however, we may suggest a smaller diameter fiber.
- What is the entrance aperture size of your spectrometer? Make sure that your fiber diameter size and the entrance aperture to your spectrometer are compatible and are configured properly for your application needs.
- 3. If you have too much light in your setup, are there ways you can attenuate the light? We believe that it's better to have too much light than not enough.

Optical Fiber Diameter Sizes Available for Custom Assemblies

- 8 μm
- 400 μm
- 50 μm
- 455 μm
- 115 μm
- 500 μm
- 200 μm
- $550~\mu m$
- 230 μm
- 600 μm
- 300 μm320 μm
- 800 μm
 1000 μm

Custom Option: Jacketing Options

The fiber assembly jacketing is designed to protect the fiber and provide strain relief. But we have jacketing options that can do so much more. We offer over 15 different jacketing options; our most popular selections are listed below. With options from PEEK polymer jacketing (No. 3) designed for environments with temperatures up to 260 °C to stainless-steel BX cable with neoprene synthetic rubber (No. 8) designed for excellent mechanical tolerance, we are sure to have the right jacketing for your application environment. Each jacketing has a set of specifications so that you and one of our Applications Scientists can choose the best option for your setup.



| Ite | m | Description | Temp. Limits | Chemical Resistance | Steam Sterilizable | Mechanical Tolerance | Length Limits |
|-----|--------------------------------------|---|-----------------|------------------------|-----------------------|-------------------------|------------------|
| 1 | PVC Monocoil | OEM applications only | 70 °C | Poor | No | Good | 6 m |
| 2 | Zip Tube Blue PVDF | Best for budget-conscious applications; standard in Laboratory-grade Assemblies | 100 °C | Poor | No | Good | 50 m |
| 3 | PEEK | Good for hydrogen peroxide applications | 260 °C | Excellent | Yes | Good | 10 m |
| 4 | Zip Tube Blue PVDF | Best for budget-conscious applications; larger diameter than #2 | 100 °C | Poor | No | Good | 50 m |
| 5 | PVC over Brass BX Tube | Glossy PVC covered over brass BX tube | 100 °C | Good | No | Good | 20 m |
| 6 | Silicone Monocoil | High-end jacketing; standard in Premium- grade Assemblies | 250 °C | Good | Yes | Good | 20 m |
| 7 | Stainless-steel BX | OEM applications only; optional polyolefin heatshrink overcoat | 250 °C | Good | Yes | Poor | 4 m |
| 8 | Stainless-steel fully interlocked BX | Excellent stainless steel jacketing supports longer lengths of fiber; optional polyolefin heatshrink overcoat | 250 °C | Good | Yes | Excellent | 40 m |
| 9 | Brake cable | OEM applications only | 120 °C | Good | No | Good | 20 m |
| 10 | Stainless-steel BX with heat shrink | Neoprene synthetic rubber over stainless steel; resists heat | 120 °C | Good | No | Excellent | 10 m |
| 11 | Large-diameter stainless-steel BX | OEM applications only; optional polyolefin heatshrink overcoat | 250 °C | Good | Yes | Poor | 4 m |

^{*} Custom printing available for OEMs.

Custom Fiber & Probe Assemblies

Custom Option: Connectors & Connector Adapters

Our fiber assemblies are available with several connector options. For an upgrade fee that includes the cost of the custom connector and labor, we will replace the standard SMA 905

Connector (included in the assembly price) with any custom connector from the list below (priced separately). When ordering custom connectors, please specify the diameter size of the optical fiber to which it will be attached. You also can order connectors separately.



| Item | Description | Connector Price |
|-------------|---|--------------------|
| CONN-ST | Stainless-steel ST Connector | \$20 |
| CONN-FC | Stainless-steel FC Connector | \$18 |
| CONN-QSMA | Premium-grade SMA 905 Connector (standard in Premium-grade assemblies) | \$16 |
| CONN-SMA | Laboratory-grade SMA 905 Connector (standard in Laboratory-grade assemblies) | \$13 |
| CONN-PSMA | Process-grade SMA 905 Connector used with an assembly with Tefzel jacketing | \$23 |
| CONN-QSMA-O | Premium-grade SMA 905 Connector with O-ring | \$26 |
| CONN-SMA-O | Laboratory-grade SMA 905 Connector with O-ring | \$23 |
| CONN-PSMA-O | Process-grade SMA 905 Connector with O-ring | \$33 |
| CONN-LSMA | Laser SMA 905 Connector for use during laser or other high-intensity applications | \$30 |

Custom Option: Connector Adapters

Connector adapters allow you to mate an item with an SMA 905 Connector to an item with either an ST or FC Connector.

| Item | Description | Price |
|------------|--|-------|
| SMA-ST-ADP | SMA-to-ST Adapter for interfacing an item with an SMA 905 Connector to | \$50 |
| | an item that has an ST Connector | |
| SMA-FC-ADP | SMA-to-FC Adapter for interfacing an item with an SMA 905 Connector to | \$50 |
| | an item that has an FC Connector | |





Custom Option: Ferrules for Probe Assemblies











| | Description |
|---|--|
| 1 | 1/4" diameter stainless-steel ferrule often used in solution transmission measurements |
| 2 | 1/4" diameter PEEK ferrule used in harsh environments for solution transmission measurements |
| 3 | 1/4" diameter stainless-steel ferrule used in reflection measurements |
| 4 | 1/4" diameter Torlon ferrule with cap |
| 5 | 1/4" diameter PEEK ferrule used in harsh environments |
| 6 | 1/8" diameter stainless-steel ferrule |
| | 1/16" diameter stainless-steel ferrule |
| | 1/4" diameter stainless-steel ferrule with the tip angled to 30° |
| | Fiber-to-lens ferrule that comes with a collimating lens |

Custom Option: Epoxy

| Item | Description | Operating Temp. Continuous | Operating Temp. Intermittent | Chemical Compatibility |
|---------------|--|----------------------------|---------------------------------|---------------------------|
| EPO-TEK 353ND | Standard epoxy in all fiber assemblies | 220 °C | 350 °C | Good |
| EPO-TEK 354ND | Slightly lower curing stress on the fiber than 353ND | 200 °C | 300 °C | Good |
| EPO-TEK OM125 | Lowest curing stress on the fiber of the three epoxies | 150 °C | 250 °C | Fair |
| EPOXY-TEST | Free sample slide of all epoxies for compatibility testing | NA | NA | NA |



Premium-grade Assemblies

Premium-grade Patch Cord Assemblies

Our Premium-grade Optical Fiber Assemblies are durable, high-quality assemblies that consistently deliver uniform results with minimal signal variance. These assemblies are available in a wide variety of off-the-shelf configurations. With every order, you receive a Quality Control Report that includes both the serial number and transmission curve of the assembly. Our 2-meter Premium-grade Patch Cord Assemblies are terminated with precision SMA 905 Connectors. They connect easily to spectrometers, light sources and sampling accessories and are available in standard lengths or can be ordered in custom lengths.



| Item Code | Wavelength Range | Core Diameter | Buffer/ Coating | Assembly Length | Jacketing | LTBR* | STBR** | Price |
|-----------------|---------------------|-----------------|--------------------|--------------------|--------------------|-------|--------|-------|
| QP8-2-SMA | 400-2500 nm | 8 μm ± 0.5 μm | acrylate | 2 meters | silicone monocoil | 4 cm | 2 cm | \$159 |
| QP50-2-UV/BX | 300-1100 nm | 50 μm ± 5 μm | polyimide | 2 meters | stainless-steel BX | 4 cm | 2 cm | \$149 |
| QP50-2-UV-VIS | 300-1100 nm | 50 μm ± 5 μm | polyimide | 2 meters | silicone monocoil | 4 cm | 2 cm | \$150 |
| QP50-2-VIS-NIR | 400-2500 nm | 50 μm ± 5 μm | polyimide | 2 meters | silicone monocoil | 4 cm | 2 cm | \$150 |
| QP100-2-UV/BX | 300-1100 nm | 100 μm ± 3 μm | polyimide | 2 meters | stainless-steel BX | 4 cm | 2 cm | \$149 |
| QP100-2-UV-VIS | 300-1100 nm | 100 μm ± 3 μm | polyimide | 2 meters | silicone monocoil | 4 cm | 2 cm | \$150 |
| QP100-2-VIS/BX | 400-2500 nm | 100 μm ± 3 μm | polyimide | 2 meters | stainless-steel BX | 4 cm | 2 cm | \$149 |
| QP100-2-VIS-NIR | 400-2500 nm | 100 μm ± 3 μm | polyimide | 2 meters | silicone monocoil | 4 cm | 2 cm | \$150 |
| QP200-2-SR/BX | 200-1100 nm | 200 μm ± 4 μm | polyimide | 2 meters | stainless-steel BX | 8 cm | 2 cm | \$179 |
| QP200-2-UV/BX | 300-1100 nm | 200 μm ± 4 μm | polyimide | 2 meters | stainless-steel BX | 8 cm | 4 cm | \$149 |
| QP200-2-UV-VIS | 300-1100 nm | 200 μm ± 4 μm | polyimide | 2 meters | silicone monocoil | 8 cm | 4 cm | \$150 |
| QP200-2-VIS/BX | 400-2500 nm | 200 μm ± 4 μm | polyimide | 2 meters | stainless-steel BX | 8 cm | 4 cm | \$149 |
| QP200-2-VIS-NIR | 400-2500 nm | 200 μm ± 4 μm | polyimide | 2 meters | silicone monocoil | 8 cm | 4 cm | \$150 |
| QP230-0.25-XSR | 180-900 nm | 230 μm ± 10 μm | aluminum | 25 centimeters | stainless-steel BX | 4 cm | 2 cm | \$149 |
| QP230-1-XSR | 180-900 nm | 230 μm ± 10 μm | aluminum | 1 meter | stainless-steel BX | 4 cm | 2 cm | \$179 |
| QP230-2-XSR | 180-900 nm | 230 μm ± 10 μm | aluminum | 2 meters | stainless-steel BX | 4 cm | 2 cm | \$199 |
| QP300-1-SR | 200-1100 nm | 300 μm ± 6 μm | polyimide | 1 meter | silicone monocoil | 12 cm | 6 cm | \$150 |
| QP400-025-SR | 200-1100 nm | 400 μm ± 8 μm | polyimide | 25 centimeters | silicone monocoil | 16 cm | 8 cm | \$119 |
| QP400-025-SR/BX | 200-1100 nm | 400 μm ± 8 μm | polyimide | 25 centimeters | stainless-steel BX | 16 cm | 8 cm | \$119 |
| QP400-1-UV-VIS | 300-1100 nm | 400 μm ± 8 μm | polyimide | 1 meter | silicone monocoil | 16 cm | 8 cm | \$150 |
| QP400-2-SR | 200-1100 nm | 400 μm ± 8 μm | polyimide | 2 meters | silicone monocoil | 16 cm | 8 cm | \$179 |
| QP400-2-SR/BX | 200-1100 nm | 400 μm ± 8 μm | polyimide | 2 meters | stainless-steel BX | 16 cm | 8 cm | \$179 |
| QP400-2-UV/BX | 300-1100 nm | 400 μm ± 8 μm | polyimide | 2 meters | stainless-steel BX | 16 cm | 8 cm | \$169 |
| QP400-2-UV-VIS | 300-1100 nm | 400 μm ± 8 μm | polyimide | 2 meters | silicone monocoil | 16 cm | 8 cm | \$169 |
| QP400-2-VIS/BX | 400-2500 nm | 400 μm ± 8 μm | polyimide | 2 meters | stainless-steel BX | 16 cm | 8 cm | \$169 |
| QP400-2-VIS-NIR | 400-2500 nm | 400 μm ± 8 μm | polyimide | 2 meters | silicone monocoil | 16 cm | 8 cm | \$169 |
| QP450-0.25-XSR | 180-900 nm | 455 μm ± 10 μm | aluminum | 25 centimeters | stainless-steel BX | 8 cm | 4 cm | \$159 |
| QP450-1-XSR | 180-900 nm | 455 μm ± 10 μm | aluminum | 1 meter | stainless-steel BX | 8 cm | 4 cm | \$199 |
| QP450-2-XSR | 180-900 nm | 455 μm ± 10 μm | aluminum | 2 meters | stainless-steel BX | 8 cm | 4 cm | \$239 |
| QP600-025-SR | 200-1100 nm | 600 μm ± 10 μm | polyimide | 25 centimeters | silicone monocoil | 24 cm | 12 cm | \$129 |
| QP600-025-SR/BX | 200-1100 nm | 600 μm ± 10 μm | polyimide | 25 centimeters | stainless-steel BX | 24 cm | 12 cm | \$119 |
| QP600-025-UV | 300-1100 nm | 600 μm ± 10 μm | polyimide | 25 centimeters | silicone monocoil | 24 cm | 12 cm | \$119 |
| QP600-025-VIS-N | 400-2500 nm | 600 μm ± 10 μm | polyimide | 25 centimeters | silicone monocoil | 24 cm | 12 cm | \$119 |
| QP600-1-SR | 200-1100 nm | 600 μm ± 10 μm | polyimide | 1 meter | silicone monocoil | 24 cm | 12 cm | \$179 |
| QP600-1-SR/BX | 200-1100 nm | 600 μm ± 10 μm | polyimide | 1 meter | stainless-steel BX | 24 cm | 12 cm | \$179 |
| QP600-1-UV-VIS | 300-1100 nm | 600 μm ± 10 μm | polyimide | 1 meter | silicone monocoil | 24 cm | 12 cm | \$179 |
| QP600-2-SR | 200-1100 nm | 600 μm ± 10 μm | polyimide | 2 meters | silicone monocoil | 24 cm | 12 cm | \$219 |
| QP600-2-SR/BX | 200-1100 nm | 600 μm ± 10 μm | polyimide | 2 meters | stainless-steel BX | 24 cm | 12 cm | \$219 |
| QP600-2-UV/BX | 300-1100 nm | 600 μm ± 10 μm | polyimide | 2 meters | stainless-steel BX | 24 cm | 12 cm | \$209 |
| QP600-2-UV-VIS | 300-1100 nm | 600 μm ± 10 μm | polyimide | 2 meters | silicone monocoil | 24 cm | 12 cm | \$209 |
| QP600-2-VIS/BX | 400-2500 nm | 600 μm ± 10 μm | polyimide | 2 meters | stainless-steel BX | 24 cm | 12 cm | \$209 |
| QP600-2-VIS-NIR | 400-2500 nm | 600 μm ± 10 μm | polyimide | 2 meters | silicone monocoil | 24 cm | 12 cm | \$209 |
| QP1000-2-UV/BX | 300-1100 nm | 1000 μm ± 20 μm | acrylate | 2 meters | stainless-steel BX | 30 cm | 15 cm | \$359 |
| QP1000-2-UV-VIS | 300-1100 nm | 1000 μm ± 20 μm | acrylate | 2 meters | silicone monocoil | 30 cm | 15 cm | \$359 |
| QP1000-2-VIS/BX | 400-2500 nm | 1000 μm ± 20 μm | acrylate | 2 meters | stainless-steel BX | 30 cm | 15 cm | \$359 |
| QP1000-2-VIS-NI | 400-2500 nm | 1000 μm ± 20 μm | acrylate | 2 meters | silicone monocoil | 30 cm | 15 cm | \$359 |


^{*} LTBR stands for Long Term Bend Radius, the bend radius allowed long term (such as for storage) before damaging the fiber.

^{**} STBR stands for Short Term Bend Radius, the momentary bend radius allowed before damaging the fiber.

Premium-grade Assemblies

Premium-grade Bifurcated Optical Fiber Assemblies

Our 2-meter Premium-grade Bifurcated Optical Fiber Assemblies are Y-shaped assemblies with two fibers of the same diameter side-by-side in the common end of the assembly. From the breakout of the assembly, the two fibers diverge into two legs, which can be UV-VIS, VIS-NIR or one of each --





| Item Code | Wavelength | Core Diameter | Buffer/ | Assembly | Jacketing | LTBR* | STBR** | Price |
|-----------------|---------------|----------------|-----------|----------|--------------------|-------|--------|-------|
| | Range | | Coating | Length | | | | |
| QBIF50-UV-VIS | 300-1100 nm | 50 μm ± 5 μm | polyimide | 2 meters | silicone monocoil | 4 cm | 2 cm | \$299 |
| QBIF200-UV-VIS | 300-1100 nm | 200 μm ± 4 μm | polyimide | 2 meters | silicone monocoil | 8 cm | 4 cm | \$299 |
| QBIF200-VIS/BX | 400-2500 nm | 200 μm ± 4 μm | polyimide | 2 meters | stainless-steel BX | 8 cm | 4 cm | \$299 |
| QBIF200-VIS-NIR | 400-2500 nm | 200 μm ± 4 μm | polyimide | 2 meters | silicone monocoil | 8 cm | 4 cm | \$299 |
| QBIF400-UV-VIS | 300-1100 nm | 400 μm ± 8 μm | polyimide | 2 meters | silicone monocoil | 16 cm | 8 cm | \$329 |
| QBIF400-VIS-NIR | 400-2500 nm | 400 μm ± 8 μm | polyimide | 2 meters | silicone monocoil | 16 cm | 8 cm | \$329 |
| QBIF600-UV-VIS | 300-1100 nm | 600 μm ± 10 μm | polyimide | 2 meters | silicone monocoil | 24 cm | 12 cm | \$369 |
| QBIF600-VIS/BX | 400-2500 nm | 600 μm ± 10 μm | polyimide | 2 meters | stainless-steel BX | 24 cm | 12 cm | \$369 |
| QBIF600-VIS-NIR | 400-2500 nm | 600 μm ± 10 μm | polyimide | 2 meters | silicone monocoil | 24 cm | 12 cm | \$369 |
| QBIF200-MIXED | 300-1100 nm & | 200 μm ± 4 μm | polyimide | 2 meters | silicone monocoil | 8 cm | 4 cm | \$299 |
| | 400-2500 nm | | | | | | | |
| QBIF400-MIXED | 300-1100 nm & | 400 μm ± 8 μm | polyimide | 2 meters | silicone monocoil | 16 cm | 8 cm | \$329 |
| | 400-2500 nm | | | | | | | |

Unjacketed Bulk Optical Fiber

We offer spooled, unjacketed optical fiber primarily for those interested in making their own assemblies. We offer various core diameters from 50 μ m to 600 μ m. To improve the fiber's strength and flexibility, we triple-coat it with a polyimide buffer before spooling. We offer all types of unjacketed optical fiber: High OH (high hydroxyl content for UV-VIS), Low OH (very low hydroxyl content for VIS-NIR) and Solarization-resistant (for UV-VIS). Each fiber type is optimized for use in a particular wavelength range.



| Item Code | Wavelength Range | Core Diameter | Buffer/ Coating | Fiber Type | Length | LTBR* | STBR** | Price / Meter |
|----------------|---------------------|------------------|--------------------|---------------|----------------------|-------|--------|------------------|
| FIBER-50-UV | 300-1100 nm | 50 μm ± 5 μm | polyimide | UV-VIS | specify meter length | 4 cm | 2 cm | \$15.78 |
| FIBER-50-VIS | 400-2500 nm | 50 μm ± 5 μm | polyimide | VIS-NIR | specify meter length | 4 cm | 2 cm | \$15.78 |
| FIBER-100-UV | 300-1100 nm | 100 μm ± 3 μm | polyimide | UV-VIS | specify meter length | 4 cm | 2 cm | \$15.78 |
| FIBER-100-VIS | 400-2500 nm | 100 μm ± 3 μm | polyimide | VIS-NIR | specify meter length | 4 cm | 2 cm | \$15.78 |
| FIBER-200-SR | 200-1100 nm | 200 μm ± 4 μm | polyimide | UV/SR-VIS | specify meter length | 4 cm | 2 cm | \$15.00 |
| FIBER-200-UV | 300-1100 nm | 200 μm ± 4 μm | polyimide | UV-VIS | specify meter length | 8 cm | 4 cm | \$12.63 |
| FIBER-200-VIS | 400-2500 nm | 200 μm ± 4 μm | polyimide | VIS-NIR | specify meter length | 8 cm | 4 cm | \$12.63 |
| FIBER-300-SR | 200-1100 nm | 300 μm ± 6 μm | polyimide | UV/SR-VIS | specify meter length | 12 cm | 6 cm | \$19.38 |
| FIBER-300-UV | 300-1100 nm | 300 μm ± 6 μm | polyimide | UV-VIS | specify meter length | 12 cm | 6 cm | \$18.75 |
| FIBER-300-VIS | 400-2500 nm | 300 μm ± 6 μm | polyimide | VIS-NIR | specify meter length | 12 cm | 6 cm | \$18.75 |
| FIBER-400-SR | 200-1100 nm | 400 μm ± 8 μm | polyimide | UV/SR-VIS | specify meter length | 16 cm | 8 cm | \$28.75 |
| FIBER-400-UV | 300-1100 nm | 400 μm ± 8 μm | polyimide | UV-VIS | specify meter length | 16 cm | 8 cm | \$23.38 |
| FIBER-400-VIS | 400-2500 nm | 400 μm ± 8 μm | polyimide | VIS-NIR | specify meter length | 16 cm | 8 cm | \$23.38 |
| FIBER-500-SR | 200-1100 nm | 500 μm ± 10 μm | polyimide | UV/SR-VIS | specify meter length | 20 cm | 10 cm | \$40.00 |
| FIBER-500-UV | 300-1100 nm | 500 μm ± 10 μm | polyimide | UV-VIS | specify meter length | 20 cm | 10 cm | \$31.25 |
| FIBER-500-VIS | 400-2500 nm | 500 μm ± 10 μm | polyimide | VIS-NIR | specify meter length | 20 cm | 10 cm | \$31.25 |
| FIBER-600-SR | 200-1100 nm | 600 μm ± 10 μm | polyimide | UV/SR-VIS | specify meter length | 24 cm | 12 cm | \$20.34 |
| FIBER-600-UV | 300-1100 nm | 600 μm ± 10 μm | polyimide | UV-VIS | specify meter length | 24 cm | 12 cm | \$36.88 |
| FIBER-600-VIS | 400-2500 nm | 600 μm ± 10 μm | polyimide | VIS-NIR | specify meter length | 24 cm | 12 cm | \$36.88 |
| FIBER-1000-UV | 300-1100 nm | 1000 μm ± 20 μm | acrylate | UV-VIS | specify meter length | 30 cm | 15 cm | \$96.25 |
| FIBER-1000-VIS | 400-2500 nm | 1000 μm ± 20 μm | acrylate | VIS-NIR | specify meter length | 30 cm | 15 cm | \$96.25 |

^{*} LTBR stands for Long Term Bend Radius, the bend radius allowed long term (such as for storage) before damaging the fiber.


^{**} STBR stands for Short Term Bend Radius, the momentary bend radius allowed before damaging the fiber.



Laboratory-grade Assemblies

Laboratory-grade Patch Cord Optical Fiber Assemblies

Our Laboratory-grade Optical Fiber Assemblies offer high quality at an affordable price. These off-the-shelf assemblies come in various lengths and wavelength ranges, and in bifurcated and splitter designs. Our standard Laboratory-grade Assemblies act as both illumination and read fibers and connect easily to our spectrometers, light sources and accessories. Solarization-resistant Assemblies (these have -SR in the item code) are for applications below 300 nm. UV radiation below 300 nm degrades transmission in standard silica fibers, resulting in solarization (increased light absorption in the UV fiber that can invalidate data). For applications below 200 nm, we recommend our Premiumarade Xtreme Solarization-resistant Fiber and Probe Assemblies (see page 146).



 $600 \mu m \pm 10 \mu m$

 $600 \mu m \pm 10 \mu m$

1000 μm ± 20 μm

 $1000 \ \mu m \pm 20 \ \mu m$

300-1100 nm

400-2500 nm

300-1100 nm

400-2500 nm

| Specifications | | | | | | |
|-----------------|---|---------------------|---|--|--|--|
| Type: | Laboratory-grade Patch Cord Optical Fiber Assemblies | Connector(s): | Laboratory-grade SMA 905 Connectors | | | |
| Fiber profile: | Step-index multimode (all diameter assemblies except for (8 µm) and | Assembly jacketing | zip tube blue PVDF | | | |
| | Step-index single mode (8 µm diameter assemblies) | Temperature range: | -20 °C to 80 °C | | | |
| Fiber core: | Pure silica core | Numerical aperture: | Multimode: 0.22 ± 0.02 (yields acceptance angle of 24.8° in air); | | | |
| Fiber cladding: | Fluorine-doped silica cladding |] | Single mode: 0.14 ± 0.02 (yields acceptance angle of 14.94° in air) | | | |

polyimide

polyimide

acrylate

acrylate

10 meters

10 meters

2 meters

2 meters

P600-10-UV-VIS

P600-10-VIS-NIR

P1000-2-UV-VIS

P1000-2-VIS-NIR

24 cm

24 cm

30 cm

30 cm

12 cm

12 cm

15 cm

15 cm

\$349

\$349

\$299

\$299

LTBR stands for Long Term Bend Radius, the bend radius allowed long term (such as for storage) before damaging the fiber.

^{**} STBR stands for Short Term Bend Radius, the momentary bend radius allowed before damaging the fiber.

Laboratory-grade Assemblies

Laboratory-grade Bifurcated and Splitter Optical Fiber Assemblies

These 2-meter Bifurcated Optical Fiber Assemblies (right, top) are Y-shaped assemblies that have two fibers of the same diameter side-by-side in the common end, or the tail of the assembly. From the nexus or breakout of the assembly, the two fibers diverge into two separate legs. You may specify that both fibers in the assembly are UV-VIS, VIS-NIR or one of each -- a "mixed" bifurcated assembly.

A splitter (right, bottom) is a 2-meter, Y-shaped assembly with a stainless steel breakout located midway from the ends of the assembly. Each splitter is made up of three separate optical fibers, all of the same diameter, and epoxied at the nexus of the Y-shaped assembly. A splitter can route light from two different sources to illuminate one sample or from one source to illuminate two samples. Splitters have lower transmission efficiency than other fiber assemblies due to their design.



| Item Code | Wavelength Range | Core Diameter | Buffer/ Coating | Assembly Length | LTBR* | STBR** | Price |
|----------------------------|----------------------|-----------------|--------------------|--------------------|-------|--------|-------|
| Bifurcated Optical Fiber | Assemblies | | | | | | |
| BIF50-UV-VIS | 300-1100 nm | 50 μm ± 5 μm | polyimide | 2 meters | 4 cm | 2 cm | \$249 |
| BIF50-VIS-NIR | 400-2500 nm | 50 μm ± 5 μm | polyimide | 2 meters | 4 cm | 2 cm | \$249 |
| BIF200-UV-VIS | 300-1100 nm | 200 μm ± 4 μm | polyimide | 2 meters | 8 cm | 4 cm | \$249 |
| BIF200-VIS-NIR | 400-2500 nm | 200 μm ± 4 μm | polyimide | 2 meters | 8 cm | 4 cm | \$249 |
| BIF400-UV-VIS | 300-1100 nm | 400 μm ± 8 μm | polyimide | 2 meters | 16 cm | 8 cm | \$279 |
| BIF400-VIS-NIR | 400-2500 nm | 400 μm ± 8 μm | polyimide | 2 meters | 16 cm | 8 cm | \$279 |
| BIF400-MIXED | 300-1100 nm & | 400 μm ± 8 μm | polyimide | 2 meters | 16 cm | 8 cm | \$279 |
| | 400-2500 nm | | | | | | |
| BIF600-UV-VIS | 300-1100 nm | 600 μm ± 10 μm | polyimide | 2 meters | 24 cm | 12 cm | \$329 |
| BIF600-VIS-NIR | 400-2500 nm | 600 μm ± 10 μm | polyimide | 2 meters | 24 cm | 12 cm | \$329 |
| Splitter Optical Fiber Ass | semblies | | | | | | |
| SPLIT200-UV-VIS | 300-1100 nm | 200 μm ± 4 μm | polyimide | 2 meters | 8 cm | 4 cm | \$499 |
| SPLIT200-VIS-NIR | 400-2500 nm | 200 μm ± 4 μm | polyimide | 2 meters | 8 cm | 4 cm | \$499 |
| SPLIT400-UV-VIS | 300-1100 nm | 400 μm ± 8 μm | polyimide | 2 meters | 16 cm | 8 cm | \$499 |
| SPLIT400-VIS-NIR | 400-2500 nm | 400 μm ± 8 μm | polyimide | 2 meters | 16 cm | 8 cm | \$499 |
| Keyed SMA Optical Fibe | er Assemblies, Round | to Keyed Linear | | | | | |
| PL100-2-UV-VIS | 300-1100 nm | 100 μm ± 3 μm | polyimide | 2 meters | 4 cm | 2 cm | \$299 |
| PL100-2-VIS-NIR | 400-2500 nm | 100 μm ± 3 μm | polyimide | 2 meters | 4 cm | 2 cm | \$299 |
| PL100-2-MIXED | 300-1100 nm & | 100 μm ± 3 μm | polyimide | 2 meters | 4 cm | 2 cm | \$299 |
| | 400-2500 nm | | | | | | |
| PL200-2-MIXED | 300-1100 nm & | 200 μm ± 4 μm | polyimide | 2 meters | 8 cm | 4 cm | \$299 |
| | 400-2500 nm | | | | | | |

LTBR stands for Long Term Bend Radius, the bend radius allowed long term (such as for storage) before damaging the fiber.

^{**} STBR stands for Short Term Bend Radius, the momentary bend radius allowed before damaging the fiber.

| Specification | 18 |
|------------------------|---|
| Type: | Laboratory Grade Bifurcated, Splitter and Keyed |
| | Optical Fiber Assemblies |
| Fiber profile: | Step-index multimode |
| Fiber core: | Pure silica core |
| Fiber cladding: | Fluorine-doped silica cladding |
| Connector(s): | Laboratory-grade SMA 905 Connectors |
| Assembly jacketing: | zip tube blue PVDF |
| Breakout: | Midway of assembly at 1 meter |
| Operating temp. range: | -20 °C to 80 °C |
| Numerical aperture: | 0.22 ± 0.02 (yields acceptance angle of 24.8° in air) |



Xtreme Solarization-resistant Assemblies

High Transparency and Durability


Xtreme Solarization-resistant Optical Fiber and Probe Assemblies for spectroscopy are manufactured using a proprietary process that provides enhanced UV transmission -- signal will transmit to 180 nm -- and remarkable resistance to UV degradation, making it ideal for deep-UV applications (<300 nm). Ocean Optics is the only spectroscopy manufacturer to offer XSR Fiber.

Solarization in Typical Fibers

UV radiation below 300 nm degrades transmission in standard silica fibers, resulting in solarization (increased light absorption in the UV fiber that can invalidate data). For applications below 300 nm, we recommend Premium-grade XSR Xtreme Solarization-resistant Fiber and Probe Assemblies, which are available for immediate delivery in standard lengths. Shorter lengths for solarization-resistant assemblies maximize UV throughput. Custom lengths are available; however, we recommend speaking with our Applications Scientists before ordering.

High Production Values

Our XSR products are part of our high-quality Premium-grade line of optical fibers and probes. XSR Fiber and Probe Assemblies are robust and durable; each assembly has an aluminum coating, a stainless-steel BX jacketing, and our high-quality SMA 905 Connectors with captive end cap that protects fiber tips against scratches and contaminants. And like all of our optical fibers, the XSR assemblies are precisely polished to work with our miniature fiber optic spectrometers and accessories.







| Specification | าธ |
|------------------------|---|
| Type: | Premium Grade Xtreme XSR Solarization-resistant |
| | Optical Fiber Assemblies and |
| | Reflection/Backscattering Probes |
| Wavelength range: | 180-900 nm |
| Fiber profile: | Step-index multimode |
| Fiber core: | Pure silica core |
| Fiber cladding: | Fluorine-doped silica cladding |
| Fiber buffer/coating: | Aluminum |
| Operating temp. range: | -50 °C to 250 °C |
| Numerical aperture: | 0.22 ± 0.02 (yields acceptance angle of 24.8° in air) |
| Connector(s): | Premium-grade SMA 905 Connectors |
| Assembly jacketing | Stainless-steel BX |

| Item Code | Product Type | Assembly Length | Core Diameter | LTBR* | STBR** | Price |
|----------------|---------------------------------|-----------------|---------------|-------|--------|-------|
| QP230-0.25-XSR | Optical Fiber Assembly | 25 centimeters | 230 μm | 4 cm | 2 cm | \$149 |
| QP230-1-XSR | Optical Fiber Assembly | 1 meter | 230 µm | 4 cm | 2 cm | \$179 |
| QP230-2-XSR | Optical Fiber Assembly | 2 meters | 230 µm | 4 cm | 2 cm | \$199 |
| QP450-0.25-XSR | Optical Fiber Assembly | 25 centimeters | 450 µm | 8 cm | 4 cm | \$159 |
| QP450-1-XSR | Optical Fiber Assembly | 1 meter | 450 μm | 8 cm | 4 cm | \$199 |
| QP450-2-XSR | Optical Fiber Assembly | 2 meters | 450 μm | 8 cm | 4 cm | \$239 |
| QR230-7-XSR/BX | Reflection/backscattering Probe | 2 meters | 230 µm | 4 cm | 2 cm | \$599 |
| QR450-7-XSR | Reflection/backscattering Probe | 2 meters | 450 μm | 8 cm | 4 cm | \$999 |

LTBR stands for Long Term Bend Radius, the bend radius allowed long term (such as for storage) before damaging the fiber.

^{**} STBR stands for Short Term Bend Radius, the momentary bend radius allowed before damaging the fiber.

Flame Loop Fiber Optic Probe

Heat-resistant Fiber Optic Probe

The FL-400 Flame Loop Fiber Optic Probe couples to our spectrometers to measure in situ emission spectra of samples such as dissolved metals and high-temperature plasmas. The FL-400 consists of a high-temperature 400 μ m gold-jacketed UV-VIS optical fiber in an 8-inch-long nickel sleeve. The assembly operates in temperatures up to 700 °C. The probe connects to the 21-02 Splice Bushing and a P400-2-UV-VIS Optical Fiber, which couples to a spectrometer to measure emission spectra. (Components are sold separately.)

Use as a Flame Loop Probe or Use as a Heat-resistant Fiber Probe

The FL-400 is especially beneficial as an emission spectroscopy teaching tool to observe atomic emission lines of dissolved metals. You simply dip the loop in your sample material and pass the loop over an open flame to take emission measurements. To use the FL-400 as a heat-resistant emission probe, remove the flame loop and insert the FL-400 into a high-temperature environment to monitor emission.

FL-400: \$499 P400-2-UV-VIS: \$119 21-02: \$13

| Probe Assembly Specifications | | | |
|-------------------------------|--|--|--|
| Fiber diameter: | 400 μm core diameter | | |
| Fiber core/cladding: | Fused silica core and doped, fused silica cladding | | |
| Fiber jacketing: | Gold | | |
| Fiber type: | 1 single-strand, multimode fiber | | |
| Wavelengths covered: | 300-1100 nm | | |
| Probe sleeve (ferrule): | Nickel | | |
| Probe dimensions: | 17.78 cm length, 20-gauge probes with 0.902 mm OD | | |
| Temperature range: | -269 °C to 700 °C | | |
| Numerical aperture: | 0.22 | | |
| Connector: | SMA 905 | | |
| | | | |






Cosine-corrected Irradiance Probe



CC-3 and CC-3-UV Cosine Correctors collect radiation from a 180° solid angle. When screwed onto the end of an optical fiber, the cosine corrector and optical fiber become an irradiance probe, measuring the intensity of light normal to the probe surface defined by the diffusing material. The probe then couples to one of our spectrometers to make a complete spectroradiometer for relative and absolute spectral intensity measurements, such as measuring UV-A and UV-B in natural solar environments, evaluating emissive color sources and analyzing light sources such as LEDs and lasers. The CC-3 has an opaline glass diffuser for VIS-NIR; the CC-3-UV utilizes Spectralon for UV-NIR. Each disc sits flush at the end of 6.35-mm outer diameter barrel, which is threaded on one side for SMA 905 Connectors.


CC-3: \$99 CC-3-UV: \$129



| Specifications | | | |
|---------------------|---------------|-------------|--|
| | CC-3 | CC-3-UV | |
| Diffusing material: | Opaline glass | Spectralon | |
| Wavelength range: | 350-1000 nm | 200-1100 nm | |
| Disc thickness: | 7.9 mm | 7.9 mm | |
| Dimensions: | 6.35 mm OD | 6.35 mm OD | |
| Field of view: | 180° | 180° | |
| Connector: | SMA 905 | SMA 905 | |

Reflection/Backscattering Probes

Our R-series Fiber Optic Reflection Probes are used for measuring specular or diffuse reflectance from a surface, fluorescence from solid surfaces, or backscattering and fluorescence in solutions and powders. These probes come in all four fiber types (XSR, UV/SR, UV-VIS and VIS-NIR), or a combination of fiber types. (See page 139 for more on each fiber type.)



Our standard reflection probes use a 6-around-1 closepacking design to ensure parallel orientation of the fibers.

Standard Reflection/Backscattering Probes

Stainless Steel Ferrule

Our Standard Reflection/Backscattering Probes are based on a natural close-packing arrangement of optical fibers. Typically, the arrangement is a tight bundle of seven optical fibers -- six illumination fibers around one read fiber. This arrangement ensures parallel orientation of the fibers. The center or read fiber splits from the other six fibers and couples to a spectrometer. The outer six illumination fibers connect to the light source.



The RP200-7-UV-VIS consists of a bundle of 200-µm fibers in a six-around-one design, but has a PEEK ferrule for applications where samples may be corrosive to the standard stainless-steel ferrule.



The R200-REF consists of an R200-7 and an additional fiber to monitor the illumination (or reference) source, which is useful for any experiment in which variation or instability in the spectral output of the light source needs to be monitored.



The R200-MIXED has 14 fibers -- six UV-VIS and six VIS-NIR illumination fibers, plus one UV-VIS and one VIS-NIR read fiber (see bundle photo at right). It couples easily to a dualchannel spectrometer in which each channel is set for a different wavelength range.

Angled Probes for Solutions & Powders

Our angled probes also have a six-around-one design, but employ a 30° window to remove specular reflection effects when the probe is immersed in dense solutions and powders. We offer the angled probe with two different diameter fibers, 200 μ m and 400 μ m. These backscattering probes are especially useful for fluorescence measurements.





Reflection/Backscattering Probes


| Item Code | Wavelength Range | Core Diameter | Fiber Bundle | Fiber Jacketing | Probe Ferrule | LTBR* | STBR** | Price |
|-------------------|------------------------------|-------------------|---|--------------------|-----------------------------|-------|--------|-------|
| Standard Reflecti | on/Backscatter | ring Probes | | | | | | |
| QR200-7-UV-VIS | 300-1100 nm | 200 μm ± 4 μm | 6 illumination fibers around 1 read | silicone monocoil | 6.35 mm OD | 8 cm | 4 cm | \$399 |
| QR200-7-VIS-NIR | 400-2500 nm | 200 μm ± 4 μm | 6 illumination fibers around 1 read | silicone monocoil | 6.35 mm OD | 8 cm | 4 cm | \$399 |
| QR230-7-XSR/BX | 180-900 nm | 230 μm | 6 illumination fibers around 1 read, with aluminum coating | stainless-steel BX | 6.35 mm OD | | | \$599 |
| QR400-7-SR | 200-1100 nm | 400 μm ± 8 μm | 6 illumination fibers around 1 read | silicone monocoil | 6.35 mm OD | 16 cm | 8 cm | \$699 |
| QR400-7-SR/BX | 200-1100 nm | 400 μm ± 8 μm | 6 illumination fibers around 1 read | stainless-steel BX | 6.35 mm OD | 16 cm | 8 cm | \$549 |
| QR400-7-UV/BX | 300-1100 nm | 400 μm ± 8 μm | 6 illumination fibers around 1 read | stainless-steel BX | 6.35 mm OD | 16 cm | 8 cm | \$549 |
| QR400-7-UV-VIS | 300-1100 nm | 400 μm ± 8 μm | 6 illumination fibers around 1 read | silicone monocoil | 6.35 mm OD | 16 cm | 8 cm | \$650 |
| QR400-7-VIS/BX | 400-2500 nm | 400 μm ± 8 μm | 6 illumination fibers around 1 read | stainless-steel BX | 6.35 mm OD | 16 cm | 8 cm | \$549 |
| QR400-7-VIS-NIR | 400-2500 nm | 400 μm ± 8 μm | 6 illumination fibers around 1 read | silicone monocoil | 6.35 mm OD | 16 cm | 8 cm | \$649 |
| QR450-7-XSR | 180-900 nm | 450 μm | 6 illumination fibers around 1 read, with aluminum coating | stainless-steel BX | 6.35 mm OD | | | \$999 |
| QR600-7-SR-125F | 200-1100 nm | 600 μm ± 10 μm | 6 illumination fibers around 1 read | silicone monocoil | 3.18 mm OD | 24 cm | 12 cm | \$899 |
| QR600-7-UV-125F | 300-1100 nm | 600 μm ± 10 μm | 6 illumination fibers around 1 read | silicone monocoil | 3.18 mm OD | 24 cm | 12 cm | \$849 |
| QR600-7-VIS-125 | 400-2500 nm | 600 μm ± 10 μm | 6 illumination fibers around 1 read | silicone monocoil | 3.18 mm OD | 24 cm | 12 cm | \$849 |
| R200-7-UV-VIS | 300-1100 nm | 200 μm ± 4 μm | 6 illumination fibers around 1 read | zip tube blue PVDF | 6.35 mm OD | 8 cm | 4 cm | \$399 |
| R200-7-VIS-NIR | 400-2500 nm | 200 μm ± 4 μm | 6 illumination fibers around 1 read | zip tube blue PVDF | 6.35 mm OD | 8 cm | 4 cm | \$399 |
| R400-7-SR | 200-1100 nm | 400 μm ± 8 μm | 6 illumination fibers around 1 read | zip tube blue PVDF | 6.35 mm OD | 16 cm | 8 cm | \$549 |
| R400-7-UV-VIS | 300-1100 nm | 400 μm ± 8 μm | 6 illumination fibers around 1 read | zip tube blue PVDF | 6.35 mm OD | 16 cm | 8 cm | \$499 |
| R400-7-VIS-NIR | 400-2500 nm | 400 μm ± 8 μm | 6 illumination fibers around 1 read | zip tube blue PVDF | 6.35 mm OD | 16 cm | 8 cm | \$499 |
| R600-7-SR-125F | 200-1100 nm | 600 μm ± 10 μm | 6 illumination fibers around 1 read | zip tube blue PVDF | 3.18 mm OD | 24 cm | 12 cm | \$699 |
| R600-7-UV-125F | 300-1100 nm | 600 μm ± 10 μm | 6 illumination fibers around 1 read | zip tube blue PVDF | 3.18 mm OD | 24 cm | 12 cm | \$699 |
| R600-7-VIS-125F | 400-2500 nm | 600 μm ± 10 μm | 6 illumination fibers around 1 read | zip tube blue PVDF | 3.18 mm OD | 24 cm | 12 cm | \$699 |
| RP200-7-UV-VIS | 300-1100 nm | 200 μm ± 4 μm | 6 illumination fibers around 1 read | zip tube blue PVDF | PEEK ferrule, 6.35 mm OD | 8 cm | 4 cm | \$499 |
| Reflection/Backs | cattering Probe | s with Reference | Leg | | | | | |
| QR200-REF-UV-VI | 300-1100 nm | 200 μm ± 4 μm | 6 illumination fibers around 1 read plus 1 fiber to monitor illumination | silicone monocoil | 6.35 mm OD | 8 cm | 4 cm | \$549 |
| QR200-REF-VIS-N | 400-2500 nm | 200 μm ± 4 μm | 6 illumination fibers around 1 read plus 1 fiber to monitor illumination | silicone monocoil | 6.35 mm OD | 8 cm | 4 cm | \$549 |
| R200-REF-UV-VIS | 300-1100 nm | 200 μm ± 4 μm | 6 illumination fibers around 1 read plus 1 fiber to monitor illumination | zip tube blue PVDF | 6.35 mm OD | 8 cm | 4 cm | \$549 |
| R200-REF-VIS-NI | 400-2500 nm | 200 μm ± 4 μm | 6 illumination fibers around 1 read plus 1 fiber to monitor illumination | zip tube blue PVDF | 6.35 mm OD | 8 cm | 4 cm | \$549 |
| Reflection/Backs | cattering Probe | es for Expanded V | Vavelength Coverage | | | | | |
| QR200-12-MIXED | 300-1100 nm & 400-2500 nm | 200 μm ± 4 μm | 6 UV-VIS & 6 VIS-NIR illumination fibers around 1 UV-VIS & 1 VIS-NIR fibers | silicone monocoil | 6.35 mm OD | 8 cm | 4 cm | \$749 |
| R200-12-MIXED | 300-1100 nm & 400-2500 nm | 200 μm ± 4 μm | 6 UV-VIS & 6 VIS-NIR illumination fibers around 1 UV-VIS & 1 VIS-NIR fibers | zip tube blue PVDF | 6.35 mm OD | 8 cm | 4 cm | \$749 |
| Angled Probes fo | r Solutions & F | Powders | | | | | | |
| QR200-ANGLE-U | 300-1100 nm | 200 μm ± 4 μm | 6 illumination fibers around 1 read | silicone monocoil | 6.35 mm OD | 8 cm | 4 cm | \$599 |
| QR200-ANGLE-V | 400-2500 nm | 200 μm ± 4 μm | 6 illumination fibers around 1 read | silicone monocoil | | 8 cm | 4 cm | \$599 |
| QR400-ANGLE-U | 300-1100 nm | 400 μm ± 8 μm | 6 illumination fibers around 1 read | silicone monocoil | 6.35 mm OD | 16 cm | 8 cm | \$749 |
| QR400-ANGLE-V | 400-2500 nm | 400 μm ± 8 μm | 6 illumination fibers around 1 read | silicone monocoil | 6.35 mm OD | 16 cm | 8 cm | \$749 |
| R200-ANGLE-UV | 300-1100 nm | 200 μm ± 4 μm | 6 illumination fibers around 1 read | zip tube blue PVDF | | | 4 cm | \$499 |
| R200-ANGLE-VIS | 400-2500 nm | 200 μm ± 4 μm | 6 illumination fibers around 1 read | zip tube blue PVDF | | | 4 cm | \$499 |
| R400-ANGLE-UV | 300-1100 nm | 400 μm ± 8 μm | 6 illumination fibers around 1 read | zip tube blue PVDF | | 16 cm | 8 cm | \$599 |
| R400-ANGLE-VIS | 400-2500 nm | 400 μm ± 8 μm | 6 illumination fibers around 1 read | - | | | 8 cm | \$599 |
| * LTPP stands for | | | | | | | | |

^{*} LTBR stands for Long Term Bend Radius, the bend radius allowed long term (such as for storage) before damaging the fiber.

^{**} STBR stands for Short Term Bend Radius, the momentary bend radius allowed before damaging the fiber.

| Specifications | | | | |
|------------------------|--|-------------------------|---|--|
| Fiber profile: | Step-index multimode | Probe ferrule material: | Stainless steel except where noted | |
| Fiber core: | Pure silica core | Probe connector(s): | Premium-grade SMA 905 Connectors for QR-series probes, | |
| Fiber cladding: | Fluorine-doped silica cladding | | Laboratory-grade SMA 905 Connectors for R-series probes | |
| Fiber buffer/coating: | Polyimide coating except where noted | Breakout: | Midway of assembly at 1 meter | |
| Operating temp. range: | -20 °C to 80 °C | Total probe length: | 2 meters for fiber assembly, 76.2 mm long for probe ferrule = | |
| Numerical aperture: | 0.22 ± 0.02 (yields an acceptance angle of 24.8° in air) | | 2.076 meters; custom lengths available | |

Transmission Dip Probes



probes are a cost-effective option for many on-line and lab applications.

In Situ Monitoring

T300-RT Design

The T300-RT-UV-VIS Transmission Dip Probe consists of two 300-µm solarization-resistant optical fibers -- one illumination optical fiber and one read optical fiber -- in a 3.175-mm (1/8") outer diameter stainless steel assembly that slides into a 127-mm long, 6.35-mm (1/4") outer diameter stainless steel ferrule. Each leg of the assembly has an SMA 905 termination so that one leg can be attached to a light source and the other to a spectrometer.

T200-RT Design

The T200-RT-VIS-NIR Transmission Dip Probe has the same optical design as the T300-RT-UV-VIS, but is made with two 200-µm VIS-NIR optical fibers in its assembly.

Process Applications

Both probes fit into a standard 1/4" Swagelok fitting for installation into a pipe or reactor. Probe optics are mounted with an epoxy that offers superior chemical resistance to most solvents and can tolerate high temperatures.

Screw-on Tips: Choose Your Pathlength

Available separately are the screw-on, interchangeable probe tips necessary to use either probe. The RT-series tips come in pathlengths of 2 mm, 4 mm, 5 mm or 10 mm so that sampling setups can be configured for optically dense or dilute solutions.

| T300-RT-UV-VIS: | \$750 |
|------------------|-------|
| T200-RT-VIS-NIR: | \$750 |
| RT-2MM: | \$240 |
| RT-4MM: | \$240 |
| RT-5MM: | \$240 |
| RT-10MM: | \$240 |



The stainless steel RT-series tips screw onto the end of the T300-RT or T200-RT.

Transmission Dip Probe Operation Second-surface mirror Illumination fiber Plano-convex lens

| Specification | 18 |
|-------------------------|---|
| Fiber type: | T300: 300 µm UV/SR fiber type (200-1100 nm) |
| | T200: 200 µm VIS-NIR fiber type (400-2500 nm) |
| Pressure limit: | 100 psi |
| Temperature limit: | 100 °C without sleeve |
| Outer diameter: | 6.35 mm |
| Probe length: | 127 mm |
| Fiber length: | 2 meters |
| Breakout: | 1.5 meters from the end of the probe |
| Optics: | Fused silica |
| Probe wetted materials: | Stainless steel, fused silica, EPO-TEK 353ND |
| Pathlength: | 2, 4, 5 or 10 mm |
| Fiber jacketing: | PVC Monocoil |
| Connector: | SMA 905 |
| Probe sleeve: | Stainless steel |

Transmission Dip Probes

Dip Probe for Hostile Environments

The TP300-UV-VIS Transmission Dip Probe couples to our spectrometers and light sources to measure the absorbance and transmission of solutions in harsh environments.

Probe Assembly

The TP300-UV-VIS consists of two $300-\mu m$ optical fibers -- one illumination optical fiber and one read optical fiber -- in a 3.175-mm (1/8") outer diameter stainless steel assembly sealed into a PEEK polymer sleeve. You have the option of choosing solarization-resistant fibers or VIS-NIR fibers for your assembly.

PEEK Polymer Sleeve

The sleeve is designed for environments with temperatures up to 200 °C. The PEEK material is also radiation-resistant and has low flammability and excellent chemical resistance. The PEEK sleeve (TPSLEEVE) comes with the TP300-UV-VIS. An additional stainless steel sleeve (T300SLEEVE) is available.


Theory of Operation

The TP300 works the same way as the T300-RT and T200-RT probes. Light travels from the light source into the illumination leg of the probe and through a lens near the end of the probe. The light then transmits through the sample compartment to a second-surface mirror. The light reflects and travels back through the sample compartment a second time and is then focused by the lens onto the read fiber and through the read leg of the probe to the spectrometer.

Adjustable-pathlength Tips

There are two adjustable-pathlength tips (2-mm to 10-mm pathlengths or 10-mm to 20-mm pathlengths) available for the TP300-UV-VIS. Additionally, an RT-PH tip for mounting pH films in the optical path can be used for pH-sensing applications.

TP300-UV-VIS: \$750
TP300-VIS-NIR: \$750
T300SLEEVE: \$250
TPSLEEVE: \$250
RTP-2-10 Tip: \$240
RTP-10-20 Tip: \$240
RT-PH Tip: \$240







The RTP-series Tips (above) are for transmission and absorbance measurements. The RT-PH Tip (left) turns the TP300 Probe into a pH Sensor (see page 76).

| Specifications | | | | | |
|-----------------|--|--------------------|---|--|--|
| Fiber type: | TP300-UV-VIS 300 µm UV/SR fiber type (200-1100 nm) | Pathlength: | Adjustable from 2-10 mm or from 10-20 mm | | |
| | TP300-VIS-NIR VIS-NIRfiber type (400-2500 nm) | | RT-PH - fixed 16-mm pathlength | | |
| Outer diameter: | 3.175 mm | Outer materials: | PVDF for jacketing, PEEK polymer for sleeve | | |
| Length: | 107.9 mm for probe, 2 meters for fiber | Temperature limit: | 200 °C with PEEK sleeve | | |
| Optics: | Fused silica | Pressure limit: | 100 psi | | |

Industrial Process Probes





| Specificati | one |
|----------------------|--|
| Fiber type: | TI300-UV-VIS 300 µm diameter |
| | UV/SR fiber type (200-1100 nm) |
| | TI300-VIS-NIR 300 µm diameter VIS-NIR |
| | fiber type (400-2500 nm) |
| Wavelength range: | TI300-UV-VIS 200-1100 nm |
| | TI300-VIS-NIR 400-2500 nm |
| Pressure limit: | 250 psi |
| Temperature limit: | 300 °C |
| Sampling tip body: | 303 stainless steel |
| Sampling tip optics: | Quartz back-coated mirror and quartz lens |
| Sampling tip O-ring: | Parker perfluoroelastomer (Parofluor ULTRA) |
| | O-ring seal |
| Probe ferrule: | 12.7 mm outer diameter 303 stainless steel |
| Probe jacketing: | Fully interlocked stainless-steel jacketing over |
| | Teflon tubing; total 6.8 mm outer diameter |
| Length: | Fiber 2 meters |
| | Ferrule 12.7 cm without tip |
| | Tips 2.6 cm to 4.99 cm, depending on tip |
| Breakout distance: | 1 meter from the end of the probe |
| Immersible length: | 12.7 cm |
| Optical pathlengths: | 2, 5, 10, 25 and 50 mm pathlengths available |
| Connectors: | SMA 905 |

Industrial Environments up to 250 psi, 300 °C

Our TI300-series Transmission Industrial Dip Probes can be used in environments with pressure limits up to 250 psi and at temperatures up to 300 °C. The TI300-UV-VIS uses 300 μ m diameter solarization-resistant optical fiber

(200-1100 nm), while the TI300-VIS-NIR uses 300 μ m diameter VIS-NIR optical fiber (400-2200 nm). The TI300 probes couple to our spectrometers and light sources to measure solutions absorbance and transmission in industrial applications.

Sampling Tips Use O-rings, Replacing Epoxy

With our other transmission probe offerings, we use high-grade epoxy to adhere the sampling optics to the

sampling tips. However, most epoxies lose their adhesive properties in continuous heat over 220 °C. With the Tl300s, we mounted the probe optics into the sampling tips using Parker perfluoroelastomer (Parofluor ULTRA) O-ring seals. The material in these special O-rings offers broad chemical resistance, excellent thermal stability and temperature resistance up to 300 °C. The other materials in the screw-on, interchangeable sampling tips are grade 303 stainless steel, a back-coated quartz mirror and a quartz lens. The tips come in pathlengths of 2, 5, 10, 25 and 50 mm so that sampling setups can be configured for optically dense or dilute solutions.

Probe Ferrule & Jacketing

The Tl300s use a fully interlocked stainless-steel jacketing over Teflon tubing and have an outer diameter of 0.68 cm. The immersible part of the probe, the ferrule, is also made from 303 stainless steel and measures 12.7 cm in length, with an outer diameter of 1.27 cm.

How it Works

In a TI300 probe, there are two 300-\$\mu\$m optical fibers -- one illumination optical fiber and one read optical fiber -- in a 12.7-mm diameter stainless-steel ferrule. Light transmits via the illumination fiber through a plano-convex lens and the sample compartment to a flat, second-surface mirror (see diagram on page 150). The light reflects from this mirror, travels back through the sample compartment and is focused by the lens onto the read fiber and through the read leg of the probe to the spectrometer. The trade-offs with these probes are that they measure both transmitted light and backscattered light from the sample and have internal reflections that limit the dynamic range of the measurement. But at less than \$1,600 (for the Tl300 probe and one sampling tip), a Tl300 probe, with its high pressure and temperature limits, is a great option for many industrial applications.

| TI300-UV-VIS: | \$1,299 |
|----------------|---------|
| TI300-VIS-NIR: | \$1,299 |
| RT-TI-2MM: | \$299 |
| RT-TI-5MM: | \$299 |
| RT-TI-25MM: | \$299 |
| RT-TI-50MM: | \$299 |

Vacuum Feedthroughs

Feedthroughs with Industry-standard Flanges

These Vacuum Feedthroughs are welded into industry standard flanges, and designed for monitoring highvacuum applications from inside a vacuum system with external equipment, such as our spectrometers. The optical fiber inside of the VFT is hermetically sealed into a stainless steel shell. These VFTs come in two standard flange types: Conflat Flanges and ISO KF Flanges. The flanges use surgical-grade stainless steel with glassceramic seals. The flange-to-seal's metal-to-glass design allows the entire assembly to perform up to either 250 $^{\circ}$ C and 10⁻¹⁰ Torr for the Conflat Flanges or 150 $^{\circ}$ C and 10⁻⁸ Torr for the ISO KF Flanges.

Range

Fiber Type



| JV-133 | ١ |
|--------|---|
| | , |
| | , |
| 1 | , |
| 7 | , |
| 1 | , |
| | , |
| | , |
| - | , |
| | , |

VFT-1000-UV-275

| 100 |
|-----|
| |
| |
| 3 |



VFT-1000-UV-16



VFT-1000-UV-40

| | VFT-200-UV-133 | 200-µm diameter UV-VIS | 300-1100 nm | 1.33" OD Conflat | \$279 |
|---|------------------|--------------------------|-------------|-------------------|-------|
| | VFT-200-VIS-133 | 200-μm diameter VIS-NIR | 400-2500 nm | 1.33" OD Conflat | \$279 |
| | VFT-400-UV-133 | 400-μm diameter UV-VIS | 300-1100 nm | 1.33" OD Conflat | \$440 |
| | VFT-400-VIS-133 | 400-μm diameter VIS-NIR | 400-2500 nm | 1.33" OD Conflat | \$440 |
| 3 | VFT-600-UV-133 | 600-µm diameter UV-VIS | 300-1100 nm | 1.33" OD Conflat | \$450 |
| | VFT-600-VIS-133 | 600-μm diameter VIS-NIR | 400-2500 nm | 1.33" OD Conflat | \$450 |
| | VFT-1000-UV-133 | 1000-µm diameter UV-VIS | 300-1100 nm | 1.33" OD Conflat | \$485 |
| | VFT-1000-VIS-133 | 1000-µm diameter VIS-NIR | 400-2500 nm | 1.33" OD Conflat | \$485 |
| | VFT-200-UV-275 | 200-μm diameter UV-VIS | 300-1100 nm | 2.73" OD Conflat | \$450 |
| | VFT-200-VIS-275 | 200-µm diameter VIS-NIR | 400-2500 nm | 2.73" OD Conflat | \$450 |
| | VFT-400-UV-275 | 400-μm diameter UV-VIS | 300-1100 nm | 2.73" OD Conflat | \$455 |
| | VFT-400-VIS-275 | 400-μm diameter VIS-NIR | 400-2500 nm | 2.73" OD Conflat | \$455 |
| | VFT-600-UV-275 | 600-μm diameter UV-VIS | 300-1100 nm | 2.73" OD Conflat | \$465 |
| | VFT-600-VIS-275 | 600-μm diameter VIS-NIR | 400-2500 nm | 2.73" OD Conflat | \$465 |
| | VFT-1000-UV-275 | 1000-µm diameter UV-VIS | 300-1100 nm | 2.73" OD Conflat | \$499 |
| 5 | VFT-1000-VIS-275 | 1000-µm diameter VIS-NIR | 400-2500 nm | 2.73" OD Conflat | \$499 |
| | VFT-200-UV-16 | 200-μm diameter UV-VIS | 300-1100 nm | 1.18" OD KF16 ISO | \$430 |
| | VFT-200-VIS-16 | 200-μm diameter VIS-NIR | 400-2500 nm | 1.18" OD KF16 ISO | \$430 |
| | VFT-400-UV-16 | 400-μm diameter UV-VIS | 300-1100 nm | 1.18" OD KF16 ISO | \$435 |
| | VFT-400-VIS-16 | 400-μm diameter VIS-NIR | 400-2500 nm | 1.18" OD KF16 ISO | \$435 |
| | VFT-600-UV-16 | 600-μm diameter UV-VIS | 300-1100 nm | 1.18" OD KF16 ISO | \$445 |
| | VFT-600-VIS-16 | 600-μm diameter VIS-NIR | 400-2500 nm | 1.18" OD KF16 ISO | \$445 |
| i | VFT-1000-UV-16 | 1000-µm diameter UV-VIS | 300-1100 nm | 1.18" OD KF16 ISO | \$480 |
| | VFT-1000-VIS-16 | 1000-μm diameter VIS-NIR | 400-2500 nm | 1.18" OD KF16 ISO | \$480 |
| | VFT-200-UV-40 | 200-μm diameter UV-VIS | 300-1100 nm | 2.16" OD KF40 ISO | \$445 |
| | VFT-200-VIS-40 | 200-μm diameter VIS-NIR | 400-2500 nm | 2.16" OD KF40 ISO | \$445 |
| | VFT-400-UV-40 | 400-μm diameter UV-VIS | 300-1100 nm | 2.16" OD KF40 ISO | \$450 |
| | VFT-400-VIS-40 | 400-μm diameter VIS-NIR | 400-2500 nm | 2.16" OD KF40 ISO | \$450 |
| | VFT-600-UV-40 | 600-μm diameter UV-VIS | 300-1100 nm | 2.16" OD KF40 ISO | \$460 |
| | VFT-600-VIS-40 | 600-μm diameter VIS-NIR | 400-2500 nm | 2.16" OD KF40 ISO | \$460 |
| | VFT-1000-UV-40 | 1000-μm diameter UV-VIS | 300-1100 nm | 2.16" OD KF40 ISO | \$495 |
|) | VFT-1000-VIS-40 | 1000-µm diameter VIS-NIR | 400-2500 nm | 2.16" OD KF40 ISO | \$495 |
| | | | | | |

VFT-series Feedthroughs

Our general-purpose VFT-series Vacuum Feedthroughs are designed to penetrate NEMA enclosures. The VFT screws into a 3/8-24 external threaded hole in the vacuum chamber, or bolts into a smooth hole with the provided nut and washer.



| Item | Fiber Type | Range | Price |
|--------------|--------------------------------|-------------|-------|
| VFT-200-SR | 200-μm diameter SR fiber | 200-1100 nm | \$299 |
| VFT-200-VIS | 200-µm diameter VIS-NIR fiber | 400-2500 nm | \$299 |
| VFT-400-SR | 400-μm diameter SR fiber | 200-1100 nm | \$299 |
| VFT-400-VIS | 400-μm diameter VIS-NIR fiber | 400-2500 nm | \$299 |
| VFT-600-SR | 600-µm diameter SR fiber | 200-1100 nm | \$299 |
| VFT-600-VIS | 600-µm diameter VIS-NIR fiber | 400-2500 nm | \$299 |
| VFT-1000-UV | 1000-µm diameter UV-VIS fiber | 300-1100 nm | \$299 |
| VFT-1000-VIS | 1000-µm diameter VIS-NIR fiber | 400-2500 nm | \$299 |

Price

| Specifications | | | |
|---------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| | VFTs with Conflat Flanges | VFTs with KF ISO Flanges | VFT-Series |
| Temperature limit: | 250 °C | 150 °C | 140 °C |
| Vacuum range: | 1 x 10 ⁻¹⁰ Torr | 1 x 10 ⁻⁸ Torr | 1 x 10 ⁻⁹ Torr |
| Numerical aperture: | 0.22, and acceptance angle of 24.8° | 0.22, and acceptance angle of 24.8° | 0.22, and acceptance angle of 24.8° |

Bare Fiber Adapter Kit



The BFA-KIT Bare Fiber Adapter Kit is for the fiber tinkerer who wants to polish bare (unjacketed) optical fiber. The kit comes with fiber polishing holders for various sizes of optical fibers.

The Bare Fiber Adapter Kit includes the following:

- 6 fiber polishing holders for various sizes of optical fiber: (1 each for 100 μ m, 200 μ m, 300 μ m, 400 μ m, $600 \, \mu \text{m}$ and $1000 \, \mu \text{m}$ optical fibers)
- a BFA-KIT-CHUCK connect-and-release adapter (which can be purchased separately as well) to fasten the SMAs onto bare optical fiber
- several pieces of wire for cleaning out the polishing holders and connect-and-release adapter

An SMA-PUCK polishing puck is not included with the BFA-KIT, but is available separately. The puck is used to polish the surface of an optical fiber.

Need a special SMA 905 Connector? We can drill out an SMA 905 Connector from 149 μ m to 2705 µm. Custom-drilled connectors are available for \$50 each.

> BFA-KIT: \$249 **BFA-KIT-CHUCK:** \$169 SMA-PUCK: \$129



The SMA-PUCK,

sold separately from the BFA-KIT.



Select the fiber polishing holder that corresponds to the diameter of your bare optical fiber.




Attach the fiber polishing holder to the front of the BFA-KIT-CHUCK connectand-release adapter.

Thread your bare optical fiber through the back of the chuck and into the holder and fasten. The holder allows you to easily work with the fiber for various purposes, such as polishing.



Fiber Tinkerer's Kit



The FT-KIT Fiber Tinkerer Kit (left) includes an assortment of randomly selected, unterminated UV-VIS and VIS-NIR optical fibers. Each fiber included in the kit will be at least one meter in length. The Fiber Termination Kit (TERM-KIT) on the following page includes all the tools needed to terminate and polish fiber.

FT-KIT: \$99

Fiber Termination Kit

Inspect & Repair Fibers

The TERM-KIT Termination Kit provides you with all the tools you need to properly polish and terminate an optical fiber. The TERM-KIT is great for inspecting, repairing and polishing optical fiber assemblies. If you would like unterminated fibers for use with the TERM-KIT, the FT-KIT Fiber Tinkerer Kit includes an assortment of optical fibers in lengths of at least one meter (see page 154 for details).

Included in Each TERM-KIT

- 4 SMA 905 Connectors for 50 μm or 100 μm fibers
- 4 SMA 905 Connectors for 200 μm optical fibers
- 4 SMA 905 Connectors for 400 μm optical fibers
- 4 SMA 905 Connectors for 600 μm optical fibers
- 4 SMA 905 Connectors for 1000 μ m optical fibers
- polishing puck
- glass polishing plate (15 cm x 15 cm)
- dozens of polishing papers
- 5-cavity crimp tool (for 2.6, 3.4, 3.8, 4.5 and 6.4 mm cavities)
- scoring tool
- inspection scope
- 2-hour cure epoxy
- optical wipes

TERM-KIT: \$599

Terminators

Our SMA 905 Connector Kits are an excellent complement to the standard Termination Kit (above). We offer six different kits, each with 10 SMA 905 Connectors of the same size, drilled for precise alignment with our optical fiber. The difference between each kit is the diameter size of the SMA 905 Connector. We offer 150 μ m, 270 μ m, 380 μ m, 490 μ m, 710 μ m and 1300 μ m diameter connectors.









| TERMKIT-QSMA-71 | 0 |
|-----------------|---|

| Item | Fiber Type | For Use with Fiber Diameters of | Price |
|-----------------|---|------------------------------------|-------|
| TERMKITSMA-150 | 10 SMA Connectors for 100-micron and 50-micron fibers | 50 μm and 100 μm | \$75 |
| TERMKITSMA-270 | 10 SMA Connectors for 200-micron fibers | 200 μm | \$75 |
| TERMKITSMA-380 | 10 SMA Connectors for 300-micron fibers | 300 μm | \$75 |
| TERMKITSMA-490 | 10 SMA Connectors for 400-micron fibers | 400 μm | \$75 |
| TERMKITSMA-710 | 10 SMA Connectors for 600-micron fibers | 600 μm | \$75 |
| TERMKITSMA-130 | 10 SMA Connectors for 1000-micron fibers | 1000 μm | \$75 |
| TERMKITQSMA-150 | 10 Premium SMA Connectors for 100-micron and 50-micron fibers | 50 μm and 100 μm | \$170 |
| TERMKITQSMA-270 | 10 Premium SMA Connectors for 200-micron fibers | 200 μm | \$170 |
| TERMKITQSMA-380 | 10 Premium SMA Connectors for 300-micron fibers | 300 μm | \$170 |
| TERMKITQSMA-490 | 10 Premium SMA Connectors for 400-micron fibers | 400 μm | \$170 |
| TERMKITQSMA-710 | 10 Premium SMA Connectors for 600-micron fibers | 600 μm | \$170 |
| TERMKITQSMA-130 | 10 Premium SMA Connectors for 1000-micron fibers | 1000 μm | \$170 |



Optical Fiber Kits

Fiber Optic Kit -- UV-VIS

We've taken our most popular laboratory-grade optical fiber assemblies and accessories and combined them into cost-saving Optical Fiber Kits -- perfect for testing, teaching or just plain tinkering. The FOP-UV Optical Fiber Kit consists of five patch cord optical fiber assemblies, the Fiber Optic Variable Attenuator, a CC-3-UV Cosine Corrector, a fiber wrench and more. By buying a kit instead of each product separately, you save over \$400. The table below lists the items in the kit.

FOP-UV KIT: \$999

| FOP-UV Items | Description | Price |
|---------------|--|-------|
| P50-2-UV-VIS | (1) 50 µm diameter optical fiber; UV-VIS | \$99 |
| P200-2-UV-VIS | (1) 200 µm diameter optical fiber; UV-VIS | \$99 |
| P400-025-SR | (2) 400 µm diameter optical fibers; UV, SR | \$198 |
| P600-2-UV-VIS | (1) 600 µm diameter optical fiber; UV-VIS | \$159 |
| 21-02 | (2) Splice bushings | \$26 |
| 21-01 | (2) Bulkhead bushings | \$18 |
| FVA-UV | (1) Fiber Optic Variable Attenuator | \$499 |
| FOT-SMAWRENCH | (1) Fiber Wrench | \$10 |
| CC-3-UV | (1) Cosine Corrector | \$129 |
| FCBARREL | (2) 6.35-mm OD barrel for FC connectors | \$58 |
| FIBER-WRAP | (3) 2-foot-long pieces of fiber wrap | \$10 |
| 74-UV | (1) Collimating Lens | \$159 |

Total, if purchased separately: \$1,46

Fiber Optic Kit -- VIS-NIR

The FOP-VIS Optical Fiber Kit consists of five patch cord optical fiber assemblies, the Fiber Optic Variable Attenuator, CC-3 Cosine Corrector, fiber wrap, a fiber wrench and more. By buying a kit instead of each product separately, you save nearly \$500. The table below lists the items included in the kit.

FOP-VIS KIT: \$999

| FOP-VIS Items | Description | Price |
|----------------|---|-------|
| P50-2-VIS-NIR | (1) 50 µm diameter optical fiber; VIS-NIR | \$99 |
| P200-2-VIS-NIR | (1) 200 µm diameter optical fiber; VIS-NIR | \$99 |
| P400-2-VIS-NIR | (2) 400 µm diameter optical fibers; VIS-NIR | \$238 |
| P600-2-VIS-NIR | (1) 600 µm diameter optical fiber; VIS-NIR | \$159 |
| 21-02 | (2) Splice bushings | \$26 |
| 21-01 | (2) Bulkhead bushings | \$18 |
| FVA-UV | (1) Fiber Optic Variable Attenuator | \$499 |
| FOT-SMAWRENCH | (1) Fiber Wrench | \$10 |
| CC-3 | (1) Cosine Corrector | \$99 |
| FCBARREL | (2) 6.35-mm OD barrel for FC connectors | \$58 |
| FIBER-WRAP | (3) 2-foot-long pieces of fiber wrap | \$10 |
| 74-UV | (1) Collimating Lens | \$159 |

Total, if purchased separately: \$1,474

The FVA-UV Fiber Optic Variable Attenuator helps control the amount of light transmitted between two fibers.




The 74-UV Collimating Lens can be screwed onto the end of a fiber to collimate light.



The FC-BARREL accepts FCterminated fibers for use with our lenses, lamps and other fixtures.

The FIBER-WRAP is used to bundle and protect optical fibers.



The 21-01 SMA Bulkhead Bushing allows you to position an optical fiber on a throughpanel such as a chamber wall.



The 21-02 Splice Bushing is an in-line adapter that mates two SMA 905 Connectors, such as those on optical fibers.

The FOT-SMAWRENCH slips over SMA 905 Connectors to secure fibers to spectrometers, collimating lenses and accessories.



Fiber & Probe Fixtures & Holders

C-Mounts

Our C-MOUNT-MIC Adapter Assembly with adjustable focusing barrel has an SMA 905 Connector in its center for attaching to optical fibers. The internal C-mount threads of this assembly allow you to adapt fiber optic spectrometers to other optical devices such as microscopes and telescopes.

The MFA-C-MOUNT also connects to optical devices such as microscopes and telescopes, but its center connector is designed to accept probes with 6.35-mm (1/4") outer diameter ferrules.

C-MOUNT-MIC: \$125 MFA-C-MOUNT: \$458







The C-MOUNT-MIC Adapter Assembly.

Phototubus Microscope Adapter

The MFA-PT Phototubus Microscope Adapter adapts to a Phototubus outlet on microscopes and accepts SMA 905-terminated optical fibers.

MFA-PT: \$424



The MFA-PT Phototubus Microscope Adapter.

Right-angle Collimating Lens Holder

The 74-90-UV is an assembly for mounting lenses at right angles, and is especially useful for applications involving awkward optical fiber routing. It has a mirror located under its cap bonded with high-temperature epoxy, and reflects light from the collimating lens to 90°. Two ports accommodate 74-series Collimating Lenses (not included).

74-90-UV: \$139



Reflection Probe Holders

The RPH-1 (far right) and RPH-2 (near right) are anodized aluminum platforms with holes drilled at 45° and 90° angles to the surface. The RPH-1 holds 6.35-mm (1/4") diameter probes but with the RPH-ADP -- an adapter that fits on the RPH-1 -- you can secure 3.17 mm (1/8") diameter probes as well. The RPH-2 is for use only with probes with SMA 905 Connectors. The Curved Surface Probe Holders accommodate 6.35-mm (1/4") outer diameter probes for measuring reflection of curved surfaces. The CSH (right) has a hole drilled at a 90° angle to the surface.

RPH-1: \$75 RPH-ADP: \$69 RPH-2: \$102 CSH: \$107 CSH-45: \$122





Optical Stages

The Single-Point Reflection Stage (at right) is a probe holder for reflection measurements of optical layers and other substrates up to 150 mm in diameter. The probe holder accommodates fiber optic probes and other sampling devices up to 6.35 mm in diameter.

The Stage-RTL-T is also a sampling system for analysis of substrate materials. The STAGE-RTL-T can be configured for reflection and transmission measurements. For details on both stages, see page 109.

\$TAGE: \$631 \$TAGE-RTL-T: \$2,303



Fiber & Probe Accessories

Bulkhead Bushing

The 21-01 SMA Bulkhead Bushing assembly is a device mount for optical fibers. The bulkhead bushing allows you to position an optical fiber on a through-panel such as a chamber wall. For example, to monitor a chamber, you could configure a sampling optic that consists of an SMA 905-terminated optical fiber screwed into the bulkhead bushing and mounted to a chamber window.

21-01: \$9



Splice Bushings

The 21-02 SMA Splice Bushings are in-line adapters that connect SMA 905-terminated optical fibers (or any two objects with SMA 905 terminations). A splice bushing consists of a 0.75" screw with female ends. The standard 21-02 is made of nickel-plated brass while the 21-02-SS is made of stainless steel. They are useful for coupling patch cords to fiber optic probes and other devices, or for any multiple-fiber application where coupling our standard optical fibers and accessories is preferable to creating costly and complex fiber optic assemblies.

21-02: \$13 21-02-SS: \$49



Bulkhead & Splice Bushing Combo

The 21-02-BH SMA Bulkhead Splice Bushing is an in-line adapter that connects SMA 905-terminated optical fibers through a chamber wall or panel. The 21-02-BH features an O-ring for sealing against the inside of the panel wall and a nut and lockwasher for mounting to the outside of the panel wall.

21-02-BH: \$23



FC Barrel

Our collimating lenses come standard with SMA 905 Connectors and interface to our SMA-terminated fibers. If you have FC-terminated fiber, you could remove the inner 6.35-mm OD SMA barrel and replace it with this FC Barrel to connect to our products.

FCBARREL: \$29



Finger Fiber Wrench

The FOT-SMAWRENCH is a wrench that slips over the hex nut of the SMA 905 Connector used in Laboratory-grade Optical Fibers and helps to easily attach the fiber to connectors on spectrometers, light sources, collimating lenses and many other accessories.

FOT-SMAWRENCH: \$10



Modemixer/Modestripper

The Modemixer/Modestripper is an in-line, 3-mm Suprasil rod that connects two SMA 905-terminated optical fibers to mix core modes and eliminate clad modes throughout 180-2100 nm.

ADP-SMA-SMA: \$268



Fiber Wrap

Fiber Safety Wrap is a protective spiral wrap used to bundle and protect optical fibers. The durable wrap comes in 2-foot sections of red, yellow and green.

FIBER-WRAP: \$10

